

Amazing Gracey™




PDT GRACEY™

More efficient contact!



LEADING COMPETITOR

Hard to adapt



- **Uniform blade curvature**
- **No sharp angulations**
- **No cornered heels**
- **Easier adaptation**
- **Fewer strokes required**

Our 'smart' methods of 100% made in the USA manufacturing enable us to significantly improve the Gracey design in a way that practitioners can easily see and experience for themselves. The blades have very uniform bends without sudden sharp angulations and cornered heels like many other manufacturers' Graceys. PDT's new Amazing Graceys™ are designed for easier adaptation to the tooth surface with more efficient contact! The practitioner uses fewer strokes resulting in less stress on both the patient and the professional and reduced patient chair-time.

Amazing Gracey™

Extended Reach Mini Graceys

These perio-designed mini-bladed instruments with extended terminal shanks are ideal for accessing deeper pockets and root planing.



Gracey 00 R001
R001R – also sold in rigid



Gracey 11-12 ER* Mini R006



Gracey 5-6 ER* Mini R003



Gracey 13-14 ER* Mini R007
R007R – also sold in rigid



Gracey 7-8 ER* Mini R004
R004R – also sold in rigid



Gracey 15-16 ER* Mini R008
R008R – also sold in rigid

Anterior Graceys



★ **Gracey 1-2 R016**
R016R – also sold in rigid



Gracey 3-4 R018



Gracey 5-6 R020

PDT's Amazing Graceys™ are anatomically color coded for faster recognition!

Pre-Molar Graceys



Gracey 7-8 R022
R022R – also sold in rigid



Gracey 9-10 R024

Posterior Graceys



★ **Gracey 11-12 R026**
R026R – also sold in rigid



★ **Gracey 13-14 R028**
R028R – also sold in rigid



Gracey 11-12 ER* R042
R042R – also sold in rigid



Gracey 13-14 ER* R043
R043R – also sold in rigid



Gracey 15-16 R030
R030R – also sold in rigid



Gracey 17-18 R032
R032R – also sold in rigid



Gracey 15-16 ER* R044

PDT's Amazing Graceys™ have uniform blade curvatures for easier adaption!

Quadrant Specific Graceys

These quadrant specific Gracey designs help with efficiency by reducing the number of times instruments need to be picked up and put down.




Gracey 11-14 R034
R034R – also sold in rigid



Gracey 12-13 R036
R036R – also sold in rigid



Gracey 11-14 ER* R046



Gracey 12-13 ER* R047

★ most popular

* extended reach - longer terminal shank