



PARADISE DENTAL TECHNOLOGIES


Creators of the Montana Jack® and Innovators of Dental Product Design



O'Hehir New Millennium™ Currettes



Trisha E. O'Hehir, RDH, BS



Get the Scoop!

- Unique scoop blades that cut a 270° perimeter
- Push, pull in any direction
- Three dimensional scaling and root planing
- No tissue trauma

These unique currettes have small scoop blades with 270° radiused cutting edges. It's like having the toe of a curette without the interference of the remaining unused portion of a traditional blade. This radical departure from traditional curette design allows the clinician to reach areas missed by traditional curette blades. Designed for access and effectiveness. Can be used with either a push or pull motion, in vertical, horizontal or diagonal directions. Perfect for flutings, furcations and line angles.



OH 1-2 R150




OH 9-10 R153




★ OH 17-18 R158



★ OH 1-2ER * R151



OH 13-14 R156



★ OH 19-20 R159



OH 7-8 R152



OH 15-16 R157



★ H5 - OH 2 R160

★ most popular

* extended reach - longer terminal shank

Push, Pull in any Direction! 3-Dimensional Scaling & Root Planing

