



PARADISE DENTAL TECHNOLOGIES

Montana Jack™ Sickle Scaler



The Montana Jack™ Sickle Scaler is 100% made in the USA by Paradise Dental Technologies. In a recent independent survey, hygienists voted the Montana Jack™ as the number one instrument they could not live without!

The Montana Jack™ has the perfect contra angle for the most universal adaptation of all instrument designs! In the anterior, the instrument adapts to provide a more favorable handle position allowing greater visibility and is not as close to the patient's face as compared to adapting an H6-7. Practitioners also report that the Montana Jack™ accesses throughout the posterior with less ergonomic stress and is more comfortable for the patient. The ease of adapting this PDT designed instrument is one of many reasons why hygienists chose the Montana Jack™ as the number one instrument over the anterior H6-7 and the posterior 204S.

The Montana Jack's™ sharp and durable blade provides convenient interproximal access with its very thin sickle blade. PDT's unique process of heat treat and cryogenics produces a molecular grain structure in the steel that results in optimal hardness without sacrificing strength. This process allows for thinner blades while maintaining the greatest durability.

The Montana Jack™ features a solid resin handle that allows for the ultimate tactile sensitivity. Weighing in at only 13 grams it is one of the lightest instruments on the market. The comfortable large diameter handle has extended knurling for superior control and torque.

The combination of these features has made the Montana Jack™ the number one instrument hygienists cannot live without. Dentists also like this instrument as an interproximal carver especially for overhangs!

Do you know Jack?

PDT, Inc.
8275 Hwy 10 W • P.O. Box 17980 • Missoula, MT 59808
LOCAL 406-626-4153 • TOLL FREE 800-240-9895 • FAX 406-626-4550
www.pdtdental.com • Look for Paradise Dental on Facebook