

Sshhh... *Finely tuned instrument at work!*




PDT's Scorpion Piezo ultrasonic scaler inserts with superior metallurgy give you quieter yet more resonant tips and provide greater cavitation. Harder steel with titanium nitride coating lasts up to four times longer. Precision hole design is inspected 100% to give you exacting mist.




- Built in visual wear indicator

Very Thin



Easy interproximal and subgingival access

Very Powerful



Strong cavitation with 40% less power

Available in nine configurations including both supra and perio tips to fit all major Piezo ultrasonic scaler units.



Made in France



PARADISE DENTAL TECHNOLOGIES

Exclusive Distributor in USA and Canada

8275 Hwy. 10 West • PO Box 17980 • Missoula, MT 59808 • 406-626-4153
Fax 406-626-4550 • 800-240-9895 • www.pdtdental.com

SCORPION PIEZO INSERTS

Supragingival Inserts

Compatible with:

T151 BSA



SATELEC, Hu-Friedy Symmetry, Vista, NSK

T161 BEM



EMS, Parkell

T171 BKP



KaVo SONOsoft, KaVo PIEZOlux

T181 BSI



SIRONA Sirosonic

T191 AIR



KaVo SONICflex

T194 AIRw



W&H, NSK Air

Perio Inserts

T153 EPSA



SATELEC, Hu-Friedy Symmetry, Vista, NSK

T163 EPEM



EMS, Parkell

T173 EPKP



KaVo SONOsoft, KaVo PIEZOlux

T183 EPSI



SIRONA Sirosonic

Insert Wrenches



- T159 CLEFBSA Fits BSA, EPSA
- T169 CLEFBEM Fits BEM, EPEM
- T179 CLEFKP Fits BKP, EPKP, AIR, AIRw
- T189 CLEFSI Fits BSI, EPSI

The dynamometric insert wrenches are torque limiting, providing safety and preventing inserts from being over tightened.



PARADISE DENTAL TECHNOLOGIES

Exclusive Distributor in USA and Canada

8275 Hwy. 10 West • PO Box 17980 • Missoula, MT 59808 • 406-626-4153
Fax 406-626-4550 • 800-240-9895 • www.pdtdental.com