



# PARADISE DENTAL TECHNOLOGIES

## Wingrove™ Titanium Implant Instruments



*Blades are solid USA medical grade titanium.*

### **R661 Wingrove™ Ti B 5-6**

This posterior solid titanium curette is designed after the standard Barnhart 5-6 contour with the more adaptable PDT blade length, radius and rounder tip specifically for adaptability, performance, and safety around **wider base** posterior implant restorations.

### **R680 Wingrove™ Ti L 3-4**

This posterior solid titanium curette is designed after the Langer 3-4 contour with the more adaptable PDT blade length, radius and rounder tip specifically for adaptability, performance, and safety around **narrow base** posterior implant restorations. Also to adapt and meet the challenges of the high water bridge and full arch cement retained implants.

### **R693 Wingrove™ Ti N 128 blunt - L 5 mini**

This anterior solid titanium curette is designed after the Nebraska 128 and the Langer 5 mini with blunted rounder tips for adaptability, performance, and safety. This instrument conforms to all the anatomy challenges in peri-implant maintenance to safely remove plaque and calculus around Hader clip bar, "O" ring ball, and Locator abutments. Also safely removes plaque and calculus from thread exposed, high water bridge, and full arch cemented anterior and posterior implants.



### *Wingrove™ "Go-To" Set*

- **T005 PDT-5 Improved Exam Kit**
- **T212 EasyView Probe 3-6-9-12**
  - or T214 3-5-7-10
  - or T216 UNC 12 (1-2-3 ~ 12)
- **R661 Wingrove™ Ti B 5-6**
- **R680 Wingrove™ Ti L 3-4**
- **R693 Wingrove™ Ti N 128 blunt - L 5 mini**

PDT, Inc.

8275 Hwy 10 W • P.O. Box 17980 • Missoula, MT 59808

LOCAL 406-626-4153 • TOLL FREE 800-240-9895 • FAX 406-626-4550

[www.pdtdental.com](http://www.pdtdental.com) • Look for Paradise Dental on FaceBook